

Manual: How to use Steadystateconc

S1. Instructions for Using ‘Steadystateconc’

Steadystateconc is a software program for calculating true steady-state metabolite concentrations by the Runge-Kutta method and the S-system method. It may be downloaded for free from <https://github.com/BioprocessdesignLab/Steadystateconc.git>. The software is programmed in Microsoft Visual Studio 2015 (**Visual C++**). However, it should work with any other C++ language software. The instructions for using the software program are as follows.

- 1) Open the program folder. It contains the main execution file (**steadystateconc.sln**) and one folder (**steadystateconc**).
- 2) Open the **steadystateconc** folder. It contains another main execution file (**steadystateconc.cpp**) and two folders (**InputFile** and **OutputFile**).
- 3) Open the **InputFile** folder to define the necessary settings for the sensitivity calculations. It contains four files (**initial.dat**, **equation.h**, **rungekutta.dat**, and **print_format.dat**). They store the following information.

InputFileFolder including the data files necessary to execute COSMOS calculations.

- equation.h	System description
- initial.dat	Data file defining the numbers of dependent and
	independent variables, the numbers of fluxes, and
	initial guesses for dependent variables and independent
	variable values, and stoichiometric matrix values.
- integration.dat	File for integrating a differential equation model with
	the Runge-Kutta method to obtain a better initial guess

- 4) Open the file **initial.dat** and input necessary values. Then save the data and close the file. As an example, consider the settings for a simple pathway metabolic reaction model shown in **Fig. S1**. **Fig. S2** displays the settings for this model.

- 5) Open the file **equation.h** and describe your system in the space bounded by “//*****”. Statements outside this range must not be changed.

- 4.1. Define variables and set their values in **double** type or **complex<double>** type in C++. Essentially any reasonable variable name can be used if it is defined here.
- 4.2. Describe the flux formulae that include the dependent variables, $x[1]$ to $x[n]$, and the independent variables, $x[n+1]$ to $x[n+m]$, in **complex<double>** C++ type.

4.3. Assign each individual flux i using the specified symbol “sv[i]=”.
Save the settings and close the file. The settings are exemplified in Fig. S3 for the simple model.




Fig. S1: Simple pathway metabolic reaction model

```

1          // //(Flag for the execution of integration using numerical methods. 0; not
           // executed)

3 1 5      // (Total number of dependent variables  $n$ ) (Total number of independent
           // variables  $m$ ) (Total number of fluxes  $l$ )
10.0       // Value of first dependent variable  $x[1]$ 
10.0       //           :           :
10.0       // Value of last dependent variable  $x[n]$ 

3.0        // Value of first independent variable  $x[n+1]$ 
           //           :           :
           // Value of last dependent variable  $x[n+m]$ 

           // Value of stoichiometric matrix  $N(X,v) \Rightarrow$ 
1   -1    0   -1    0   //  $N(x[1], sv[1]) \dots N(x[1], sv[l])$ 
0    1   -1    0    0   //           :           :
0    0    0    1   -1   //  $N(x[n], sv[1]) \dots N(x[n], sv[l])$ 

```

Fig. S2: Settings in the file **initial.dat** for the simple model shown in Fig. S1.




Fig. S3: Settings in the file **equation.h** for the simple model shown in Fig. S1.

6) The data file **integration.dat** is set for the Runge-Kutta method. The integration is executed when the first value in the data file **initial.dat** is set at 1. The format of the data is as follows.

```

0.0          // Start time
1000.0       // End of integration
0.1          // Sampling (storing) interval for calculated values ( $D=0.1$ )
100.0        // Value for devising the sampling time in order to increase the accuracy
              // of calculated values ( $d=0.1/100.0=0.001$ ).

```

7) Close the **InputFile** folder.

8) Open and run the main execution file **steadystateconc.sln**. The results are obtained very quickly.

9) Open the **OutputFile** folder to inspect the calculated results. It contains the following five files. All calculated values should be presented within round-off error.

OutputFile Folder including the files generated as a result of the calculation.

- | **RK_Conc.dat** Calculated data file, including pseudo steady-state concentrations calculated by Runge-Kutta method.
- |
- | **Steady-state_Conc.dat** Calculated data file, including true steady-state metabolite concentrations.
- |
- | **timecourse.dat** Data file for Runge-Kutta method.